


SPECIALTY INTUBATION

RUSCH® AGT ORAL PREFORMED CUFFED ETT

Description:

The Rusch® AGT Oral Preformed Cuffed Endotracheal Tube (ETT) features:

- Incorporates a preformed curvature designed to reduce tube kinking and improve surgical access by directing the tube toward the patient's chin
- Thermosensitive Ruschelit™ PVC adapts to the anatomical situation at body temperature
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
111780035	Rusch® AGT Oral Preformed Cuffed ETT, Size 3.5	10
111780040	Rusch® AGT Oral Preformed Cuffed ETT, Size 4.0	10
111780045	Rusch® AGT Oral Preformed Cuffed ETT, Size 4.5	10
111780050	Rusch® AGT Oral Preformed Cuffed ETT, Size 5.0	10
111780055	Rusch® AGT Oral Preformed Cuffed ETT, Size 5.5	10
111780060	Rusch® AGT Oral Preformed Cuffed ETT, Size 6.0	10

Catalog #	Size	Qty/Box
111780065	Rusch® AGT Oral Preformed Cuffed ETT, Size 6.5	10
111780070	Rusch® AGT Oral Preformed Cuffed ETT, Size 7.0	10
111780075	Rusch® AGT Oral Preformed Cuffed ETT, Size 7.5	10
111780080	Rusch® AGT Oral Preformed Cuffed ETT, Size 8.0	10
111780085	Rusch® AGT Oral Preformed Cuffed ETT, Size 8.5	10
111780090	Rusch® AGT Oral Preformed Cuffed ETT, Size 9.0	10

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.