

GUEDEL ORAL AIRWAYS

RUSCH® COLOR-CODED GUEDEL AIRWAY

Description:

The Rusch® Color-Coded Guedel Airway features:

- Color-coded for easy identification
- Smooth finish and rounded edges
- Reinforced bite block helps prevent occlusion
- Full range of sizes and convenient packaging options
- Emergency kit available with a clear plastic case



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
122350	0	50 mm	5	Blue	10 (Bulk)
122360	1	60 mm	6	Black	10 (Bulk)
122370	2	70 mm	7	White	10 (Bulk)
122380	3	80 mm	8	Green	10 (Bulk)
122390	4	90 mm	9	Yellow	10 (Bulk)

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
1223100	5	100 mm	10	Red	10 (Bulk)
122550	0	50 mm	5	Blue	50 (Bulk)
122560	1	60 mm	6	Black	50 (Bulk)
122570	2	70 mm	7	White	50 (Bulk)
122580	3	80 mm	8	Green	50 (Bulk)
122590	4	90 mm	9	Yellow	50 (Bulk)
1225100	5	100 mm	10	Red	50 (Bulk)
122450	0	50 mm	5	Blue	10 (Individually Pkg)
122460	1	60 mm	6	Black	10 (Individually Pkg)
122470	2	70 mm	7	White	10 (Individually Pkg)
122480	3	80 mm	8	Green	10 (Individually Pkg)
122490	4	90 mm	9	Yellow	10 (Individually Pkg)
1224100	5	100 mm	10	Red	10 (Individually Pkg)
122650	0	50 mm	5	Blue	50 (Individually Pkg)
122660	1	60 mm	6	Black	50 (Individually Pkg)

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
122670	2	70 mm	7	White	50 (Individually Pkg)
122680	3	80 mm	8	Green	50 (Individually Pkg)
122690	4	90 mm	9	Yellow	50 (Individually Pkg)
1226100	5	100 mm	10	Red	50 (Individually Pkg)
396278	0-5 (6 Airways)	50 mm-100 mm (6 Airways)	5-10 (6 Airways)	n/a	1 (Bulk)

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.