

PERIPHERAL NERVE STIMULATORS

STIMPOD NERVE STIMULATOR

Description:

The Stimpod NMS410 Peripheral Nerve Stimulator provides clinicians audiovisual feedback to facilitate accurate needle and catheter placement, and features:

- Clear, concise diagnostic screen
- Real time feedback of the actual waveform delivered
- Auto-sensing cables automatically switch stimulation current from mapping probe to needle
- Proximity Indicator provides audible and visual feedback when selected target current and pulse width ranges are reached
- Choice between linear and customizable non-linear progression of the current and pulse width settings
- Backlit screen for low light settings
- Combined nerve mapping and nerve locating cable



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Stimpod Stimulator & Accessories

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.