

Home // Pain Management // Peripheral Nerve Block Needles // StimuQuik® PNB Needle (stimulating)

PERIPHERAL NERVE BLOCK NEEDLES

STIMUQUIK® PNB NEEDLE (STIMULATING)

Description:

The clinician-inspired Arrow® StimuQuik® Stimulating Peripheral Nerve Block Needle is designed to help increase overall block success by enhancing needle control and maneuverability during needle positioning.

- Designed with less traumatic B-bevel tip for optimal balance between tactile feedback and penetration
- Insulated needle shaft with a stimulating tip allows for needle guidance using stimulation and/or ultrasound
- Centimeter markings on needle shaft for visual confirmation of insertion depth
- Ergonomic needle hub for precision grip
- Extended length of injection tubing and electrical lead connections reduces the risk of contamination to the sterile field



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

StimuQuik® PNB Needle (stimulating)