

Home // Pain Management // Peripheral Nerve Block Needles // StimuQuik® PNB Needle (stimulating)

PERIPHERAL NERVE BLOCK NEEDLES

STIMUQUIK® PNB NEEDLE (STIMULATING)

Description:

The clinician-inspired Arrow® StimuQuik® Stimulating Peripheral Nerve Block Needle is designed to help increase overall block success by enhancing needle control and maneuverability during needle positioning.

- Designed with less traumatic B-bevel tip for optimal balance between tactile feedback and penetration
- Insulated needle shaft with a stimulating tip allows for needle guidance using stimulation and/or ultrasound
- Centimeter markings on needle shaft for visual confirmation of insertion depth
- Ergonomic needle hub for precision grip
- Extended length of injection tubing and electrical lead connections reduces the risk of contamination to the sterile field



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

StimuQuik® PNB Needle (stimulating)

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.