

SPINALS

SUREBLOCK® SPINAL KITS AND SETS

Description:

The Arrow® SureBlock® Spinal Anesthesia line features the Sprotte® Spinal Needle with a variety of procedural components in kit and set configurations. With its 'atraumatic' pencil-point tip, the Sprotte® Spinal Needle has shown to have lower incidence of postdural puncture headaches as compared to the cutting-point tip needle.<sup>(1-2)</sup>

- Kits include procedural drugs and a variety of procedural and prep components
- Safety Kits include 'Kit' components with these upgrades: Safety Injection Needles and SharpsAway® II Locking Disposal Cup
- Sets include 'Kit' components, but no drugs or drug identification labels



United States

To order this product, please call:  
**866-246-6990**

Canada

To order this product, please call:  
**800-387-9699**

Lat  
Am

Kits

Safety Kits

Sets

References: 1. Vallejo MC, et al. Postdural Puncture Headache: A Randomized Comparison of Five Spinal Needles in Obstetric Patients. Anesthesia & Analgesia 2000; 91: 916–20. 2. Tung CE, et al. Cost Comparison Between the Atraumatic and Cutting Lumbar Puncture Needles. Neurology 2012; 78: 109–113.