

SPINALS

SUREBLOCK® SPINAL KITS AND SETS

Description:

The Arrow® SureBlock® Spinal Anesthesia line features the Sprotte® Spinal Needle with a variety of procedural components in kit and set configurations. With its 'atraumatic' pencil-point tip, the Sprotte® Spinal Needle has shown to have lower incidence of postdural puncture headaches as compared to the cutting-point tip needle.⁽¹⁻²⁾

- Kits include procedural drugs and a variety of procedural and prep components
- Safety Kits include 'Kit' components with these upgrades: Safety Injection Needles and SharpsAway® II Locking Disposal Cup
- Sets include 'Kit' components, but no drugs or drug identification labels



United States

To order this product, please call:
866-246-6990

Canada

To order this product, please call:
800-387-9699

Lat
Am

Kits

Safety Kits

Sets

References: 1. Vallejo MC, et al. Postdural Puncture Headache: A Randomized Comparison of Five Spinal Needles in Obstetric Patients. Anesthesia & Analgesia 2000; 91: 916–20. 2. Tung CE, et al. Cost Comparison Between the Atraumatic and Cutting Lumbar Puncture Needles. Neurology 2012; 78: 109–113.

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.