

EPIDURALS

THERACATH® EPIDURAL KITS AND SETS

Description:

The Arrow® TheraCath® Epidural Catheter is designed for pain management procedures where increased steerability is an important factor, and features a wire-reinforced, styletted body with a soft, spring-wire tip.

- Kits include procedural drugs and a variety of procedural and prep components
- Intermediate Sets include 'Kit' contents, but no drugs or drug identification labels (Note: some of these trays contain povidone-iodine prep solution)



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catheter

Intermediate Sets

Kits