

[Home](#) // [Airway Management](#) // [Humidification](#) // [Heat & Moisture Exchangers \(HME\)](#) // [Aqua+® HCH](#)

HEAT & MOISTURE EXCHANGERS (HME)

AQUA+® HCH

Description:

The Aqua+® Hygroscopic Condenser Humidifier (HCH) line uses a fiber-free, cellular, hygroscopically treated media and is available in a variety of configurations to meet a wide range of clinical needs.

- Aqua+® T and TS models, designed for spontaneously breathing tracheostomy patients, include built-in supplemental oxygen port
- Aqua+® TS includes port for convenient suctioning
- Aqua+® Flex utilizes a 1-piece construction with flexible tubing to aid in patient positioning
- Aqua+® N is ideal for neonate and infant patients
- Aqua+® 1HS and Aqua+® are designed for most adult and pediatric applications



To order this product, please call:

+1 (919) 433-4999

Catalog #	Description	Qty/Case
1570	Aqua+ Flex HCH	20
1571	Aqua+ HCH	30
1572	Aqua+ N HCH	50
1573	Aqua+ T HCH	30
1574	Aqua+ TS HCH	30

Catalog #	Description	Qty/Case
1589	Aqua+ 1HS HCH	30

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.