


**ENDOBRONCHIAL TUBES**

**DR. WHITE - RIGHT BRONCHUS ENDOBRONCHIAL TUBE**

**Description:**

The Rusch® Bronchopart™ Dr. White - Right Bronchus Endobronchial Tube features:

- Ruscheliit™ PVC technology provides the right degree of stiffness, flexibility and kink resistance
- Specially shaped bronchial cuff helps prevent accidental right upper lobe occlusion
- High-volume/low-pressure cuffs
- Color-coded for easy port identification, balloons and 15 mm connectors
- Double swivel connections with double caps
- Anatomically-correct curvatures facilitate proper placement
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Size	Qty/Box
116201350	Rusch® Dr. White - Right Bronchus Endobronchial Tube, 35 Fr	1
116201370	Rusch® Dr. White - Right Bronchus Endobronchial Tube, 37 Fr	1
116201390	Rusch® Dr. White - Right Bronchus Endobronchial Tube, 39 Fr	1
116201410	Rusch® Dr. White - Right Bronchus Endobronchial Tube, 41 Fr	1