

Home // Airway Management // Endotracheal Tubes // Endobronchial Blockers

## ENDOBONCHIAL BLOCKERS

Teleflex is proud to offer the EZ-Blocker™ Endobronchial Blocker: an intuitive, secure solution for clinicians to achieve one-lung ventilation and reduce the risk of intraoperative malpositioning. It was designed by anesthesiologists who challenged accepted clinical practices and wanted a better way to treat patients.



### RUSCH® EZ-BLOCKER™ ENDOBONCHIAL BLOCKER

The unique bifurcated distal end of the Rusch® EZ-Blocker™ Endobronchial Blocker allows for the intuitive placement of the cuffs in the right and left bronchi. Once the Rusch® EZ-Blocker™ Endobronchial Blocker is advanced through the distal end of a single-lumen endotracheal tube, the bifurcated cuffs separate and are naturally directed into the left and right main stem bronchi.

[View Details](#)

© 2023 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.