

HEAT & MOISTURE EXCHANGERS (HME)

FLEX TUBES

Description:

Teleflex offers Flex Tubes to meet a variety of clinical needs in ICU and anesthesia settings.

- Smooth-Flo® Flex Tube is supplied sterile and is ideal for head and neck surgery
- Expandi-Flex™ Flex Tube is a lightweight, positionable flex tube with rigid connectors and djustable length
- Softflex™ Flex Tube is an extremely light and flexible tube
- Individually packaged



To order this product, please call:
+1 (919) 433-4999

Catalog #	Description	Qty/Case
1800402	Softflex	100
20011	Smooth-Flo Flex	20
22362	Expandi Flex	25