

HEAT & MOISTURE EXCHANGERS (HME)

GIBECK® HUMID-FLO® HME

Description:

The exclusive design of the Gibeck® Humid-Flo® HME allows it to remain in-line during the first 72 hours of mechanical ventilation, even during aerosol treatments.

- Reduces circuit breaks aiding clinician adherence to clinical guidelines regarding infection control practices
- Protects caregivers from exposure to depressurizing circuit spray
- Allows for the delivery of SVN or MDI aerosols
- Reduces the need to routinely interrupt ventilation, therefore maintains positive-end expiratory pressure
- Hygroscopic bacteriostatic microwell paper media provide optimal moisture output with low resistance to flow
- Individually packaged



To order this product, please call:
+1 (919) 433-4999

Catalog #	Description	Qty/Case
19912	Gibeck Humid-Flo HME	25