

[Home](#) // [Airway Management](#) // [Humidification](#) // [HEPA & Depth Filters](#) // [Gibeck® ISO-Gard® Filters](#)

HEPA & DEPTH FILTERS

GIBECK® ISO-GARD® FILTERS

Description:

Gibeck® ISO-Gard® HEPA and Depth Filter Systems offer bidirectional bacterial/viral filtration for cross-contamination protection of patient, equipment and caregiver and are suitable for use on pediatric to adult patients.

- ISO-Gard® HEPA Light is an advanced, highly hydrophobic bacterial/viral filter, HEPA class 13.
- ISO-Gard® HEPA Small uses the same, proven bidirectional bacterial/viral filtration media used in the larger HEPA Light models, at half the size.
- ISO-Gard® Depth Filter is a hydrophobic, electrostatic bacterial/viral filter for use on the patient end of circuit when an HME is not required.
- Individually packaged



To order this product, please call:

+1 (919) 433-4999

Catalog #	Description	Qty/Case
19212	Depth Filter, straight	25
28002	HEPA Light with port	20
28012	HEPA Light	20
28022	HEPA Light Machine	20
28052	HEPA Small, angled	25

Catalog #	Description	Qty/Case
28062	HEPA Small, straight	25

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.