

HEAT & MOISTURE EXCHANGERS (HME)

GIBECK® TRACH-VENT® HME

Description:

The Gibeck® Trach-Vent® HME and accessory products are designed specifically for tracheostomized or intubated, spontaneously breathing patients.

- Recommended for tidal volumes of 50 mL and above
- For use on both pediatric and adult patients
- Oxy-Vent® Adaptor snaps over Trach-Vent® HME to allow delivery of supplemental oxygen
- Individually packaged



To order this product, please call:  
**+1 (919) 433-4999**

| Catalog # | Description          | Qty/Case |
|-----------|----------------------|----------|
| 41112     | Trach-Vent HME       | 50       |
| 41131     | Trach-Vent T HME     | 50       |
| 41312     | Trach-Vent+ HME      | 50       |
| 45512     | Oxy-Vent with Tubing | 10       |