

BACTERIAL/VIRAL FILTERS

MAIN FLOW BACTERIAL/VIRAL FILTER

Description:

The Main Flow Bacterial/Viral Filter features:

- 22 mm ID x 15 mm ID/22 mm OD
- Bacterial efficiency 99.999+%; viral efficiency 99.99+%
- Dead space is 45.9 mL
- Flow resistance approximately 1.8 cm H₂O at 60 LPM
- Individually packaged, 50 to a case



To order this product, please call:
+1 (919) 433-4999

Catalog #	Description	Qty/Case
1605	Main Flow Bacterial/Viral Filter, 22mm OD	50

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.