

Home // Airway Management // Laryngoscopes // MRI Laryngoscope Kit

MRI LARYNGOSCOPE KIT



RÜSCH® MRI LARYNGOSCOPE KIT

- Independently certified to support MRI conditions with a static magnetic field of 3 Tesla or less and highest gradient field of 720 Gauss/cm or less
- Unique design — black aluminum handle and gold-plated blades allow the set to be easily distinguishable
- Lightweight handle and blades
- Specially designed carbon zinc dry batteries

[View Details](#)

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.