

## RESUSCITATION

### POCKETPAK EMERGENCY MASK KIT

#### Description:

The PocketPak Resuscitator Emergency Mask Kit prevents mouth-to-mouth contact during resuscitation and fits standard manual resuscitators. The collapsible mask is easy to use and convenient to carry.

- Mask with pliable, clear dome to provide secure seal
- Capped port for delivery of supplemental oxygen
- Low-resistance valve and particle filter
- Foldable extension tube
- Packaged in a sturdy “clam shell” case, 4 masks to a convenient storage box



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.