

Home // Airway Management // LMA® Airways // Reusable LMA® Airways

REUSABLE LMA® AIRWAYS

Teleflex offers a comprehensive line of reusable LMA® Airways, providing anesthesia professionals with a full range of laryngeal mask options, so you have the right airway device for every patient.



LMA® CLASSIC® AIRWAY

- General purpose reusable airway for routine and emergency anesthetic procedures
- Silicone Cuff - Soft, flexible, and conforms to the anatomy^{1,2}
- Standard straight airway tube
- Tactile inflation pilot balloon
- Reusable up to 40 times³

[View Details](#)



LMA® PROSEAL® AIRWAY

- Versatile reusable LMA® Airway
- Double cuff design enables seal pressures of up to 32 cm H₂O¹
- Second-generation device designed to reduce risk of aspiration with an integrated drain tube and dual seal

[View Details](#)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.