

GUEDEL ORAL AIRWAYS

RUSCH® COLOR-CODED GUEDEL AIRWAY

Description:

The Rusch® Color-Coded Guedel Airway features:

- Color-coded for easy identification
- Smooth finish and rounded edges
- Reinforced bite block helps prevent occlusion
- Full range of sizes and convenient packaging options
- Emergency kit available with a clear plastic case



To order this product, please call:
+1 (919) 433-4999

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
122450	0	50 mm	5	Blue	10 (Individually Pkg)
122460	1	60 mm	6	Black	10 (Individually Pkg)
122470	2	70 mm	7	White	10 (Individually Pkg)
122480	3	80 mm	8	Green	10 (Individually Pkg)
122490	4	90 mm	9	Yellow	10 (Individually Pkg)

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
1224100	5	100 mm	10	Red	10 (Individually Pkg)
122650	0	50 mm	5	Blue	50 (Individually Pkg)
122660	1	60 mm	6	Black	50 (Individually Pkg)
122670	2	70 mm	7	White	50 (Individually Pkg)
122680	3	80 mm	8	Green	50 (Individually Pkg)
122690	4	90 mm	9	Yellow	50 (Individually Pkg)
1226100	5	100 mm	10	Red	50 (Individually Pkg)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.