

Home // Airway Management // Laryngoscopes // Single-Use Laryngoscopes // Rüsch® DispoLED® Single-Use Fiber Optic Laryngoscope Handle

## SINGLE-USE LARYNGOSCOPES

### RÜSCH® DISPOLED® SINGLE-USE FIBER OPTIC LARYNGOSCOPE HANDLE

**Description:**

LED lamp provides optimal light output

- Batteries included and easy to remove after use
- Individually packaged for single-use
- Can be tested in package without opening



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

Catalog #	Diameter	Length	Weight	Qty/Case
77700	2.8 cm	12 cm	70 g	10

© 2021 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.