

ENDOTRACHEAL TUBE ACCESSORIES

RUSCH® DOUBLE SWIVEL CONNECTOR

Description:

The Rusch® Double Swivel Connector is designed with a sealing cap to facilitate the insertion of suction catheters or fiber optic bronchoscopes, and is provided non-sterile.



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Description	Qty/Box
514900000	Rusch® Double Swivel Connector, 15 mm	25

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.