


## STYLETS

### RUSCH® SLICK STYLETTE® DEVICE

#### Description:

The Rusch® Slick Stylette® Device features:

- Made of malleable nylon-coated wire that holds its shape to facilitate intubation
- Specially formulated coating does not crack or peel, while allowing for smooth tube insertion
- Stopper helps maintain desired stylet position within endotracheal tube
- Available in three color-coded sizes for easy identification
- Single-use, sterile



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat  
Am

Catalog #	O.D. Size	Length	Qty/Box
500	6 Fr (for tube sizes 2.0-3.5)	31 cm	20
750	8 Fr (for tube sizes 4.0-6.5)	42 cm	20
1000R	10 Fr (for tube sizes 7.0-10.0)	45.5 cm	20

© 2025 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.