


STYLETS

RUSCH® SLICK STYLETTE® DEVICE

Description:

The Rusch® Slick Stylette® Device features:

- Made of malleable nylon-coated wire that holds its shape to facilitate intubation
- Specially formulated coating does not crack or peel, while allowing for smooth tube insertion
- Stopper helps maintain desired stylet position within endotracheal tube
- Available in three color-coded sizes for easy identification
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	O.D. Size	Length	Qty/Box
500	6 Fr (for tube sizes 2.0-3.5)	31 cm	20
750	8 Fr (for tube sizes 4.0-6.5)	42 cm	20
1000R	10 Fr (for tube sizes 7.0-10.0)	45.5 cm	20