


STANDARD INTUBATION

RUSCH® SUPER SAFETY™ MAGILL CUFFED ETT

Description:

The Rusch® Super Safety™ Magill Cuffed Endotracheal Tube (ETT) is a low-volume, cylindrical cuff ETT, which features:

- Ruscheliit™ PVC technology provides the right degree of stiffness, flexibility and kink resistance
- Distal tip is rounded to ease insertion
- Wide black positioning ring above cuff to help gauge depth upon insertion
- For both oral and nasal intubation
- Single-use, sterile



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Catalog #	Size	Qty/Box
112480025	Rusch® Super Safety Clear Magill Cuffed ETT, Size 2.5	10
112480030	Rusch® Super Safety Clear Magill Cuffed ETT, Size 3.0	10
112480035	Rusch® Super Safety Clear Magill Cuffed ETT, Size 3.5	10
112480040	Rusch® Super Safety Clear Magill Cuffed ETT, Size 4.0	10
112480045	Rusch® Super Safety Clear Magill Cuffed ETT, Size 4.5	10

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.