


STYLETS

SHERIDAN® SHER-I-SLIP® STYLET

Description:

The Sheridan® Sher-I-Slip® Stylet features:

- Designed to retain desired curve during oral intubation
- Material that reduces stylet-to-tube friction



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat
Am

Catalog #	O.D. Size	Length	Qty/Box
5-15101	10 Fr (for tube sizes 5.0-10.0)	34 cm	25
5-15102	5 Fr (for tube sizes 2.0-4.5)	25 cm	25
5-15103	10 Fr (for tube sizes 5.0-10.0)	51 cm	25

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.