

GUEDEL ORAL AIRWAYS

SOFT CLEAR GUEDEL AIRWAY

Description:

The Rusch® Soft Clear Guedel Airway features:

- Soft, smooth finish for maximum patient comfort
- Flexible design
- Color-coded bite block
- Full range of sizes and convenient packaging options



To order this product, please call:
+1 (919) 433-4999

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
124701002	00	40 mm	4	Pink	10 (Individually Packaged)
124701003	0	50 mm	5	Blue	10 (Individually Packaged)
124701010	1	60 mm	6	Black	10 (Individually Packaged)
124701020	2	70 mm	7	White	10 (Individually Packaged)
124701030	3	80 mm	8	Green	10 (Individually Packaged)
124701040	4	90 mm	9	Yellow	10 (Individually Packaged)

Catalog #	Size	Length	ISO Standard Size Design	color	Qty/Box
124701050	5	100 mm	10	Red	10 (Individually Packaged)
124701060	6	120 mm	12	Purple	10 (Individually Packaged)

© 2026 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.