

STANDARD INTUBATION

Teleflex's Rusch® and Sheridan® Brands of standard intubation products provides anesthesia professionals with a full range of options. We offer endotracheal tubes for short or long-term intubation, from high to low-volume, cuffed or uncuffed, so you have the right airway device for every patient.

**SHERIDAN/HVT® MURPHY EYE ETT**


The Sheridan/HVT® Murphy Eye Endotracheal Tube (ETT) is a high-volume, tapered, low pressure cuff ETT.

[View Details](#)
**SHERIDAN/CF® MURPHY EYE ETT**

The Sheridan/CF® Murphy Eye Endotracheal Tube (ETT) is a close-fitting, low-pressure cuff ETT.

[View Details](#)
**SHERIDAN/CF® MAGILL ETT**

The Sheridan/CF® Magill Endotracheal Tube (ETT) is a close-fitting, low-pressure cuff ETT.

[View Details](#)
**RUSCH® SUPER SAFETY™ MAGILL CUFFED ETT**

The Rusch® Super Safety™ Magill Cuffed Endotracheal Tube (ETT) is a low-volume, cylindrical cuff ETT.

[View Details](#)
**SHERIDAN UNCUFFED™ MURPHY EYE ETT**

The Sheridan Uncuffed™ Murphy Eye Endotracheal Tube (ETT) is a single-use ETT which softens at body temperature to conform to airway anatomy, and features:

- Smooth, molded, gently cupped tip eases insertion during intubation
- Special depth lines at distal tip reference depth of intubation

[View Details](#)
**RUSCH® FLEXI-SET™ CUFFED ETT**

The Rusch® Flexi-Set™ Cuffed Endotracheal Tube (ETT) features:

- Available with the Rusch® Safety Clear™ or the Rusch Uncuffed™ Murphy Eye



- ETT, preloaded with the Flexi-Slip™ Stylet
- Stylet has a low-friction coating to facilitate smooth withdrawal
- Preassembled and ready to use, saving time and space

[View Details](#)



RUSCH® FLEXI-SET™ UNCUFFED ETT

The Rusch® Flexi-Set™ Uncuffed Endotracheal Tube (ETT) features:

- Available with the Rusch® Safety Clear™ or the Rusch Uncuffed™ Murphy Eye ETT, preloaded with the Flexi-Slip™ Stylet
- Stylet has a low-friction coating to facilitate smooth withdrawal
- Preassembled and ready to use, saving time and space

[View Details](#)



RUSCH® SLICK SET® CUFFED ETT

The Rusch® Slick Set® Cuffed Endotracheal Tube (ETT) features:

- Available with the Rusch® Super Safety™ Clear or Rusch Uncuffed™ Murphy Eye ETT, preloaded with Slick Stylette® Device
- Preassembled and ready to use, saving time and space
- Cork-style stopper helps maintain stylet position

[View Details](#)



RUSCH® SLICK SET® UNCUFFED ETT

The Rusch® Slick Set® Uncuffed Endotracheal Tube (ETT) features:

- Available with the Rusch® Super Safety™ Clear or Rusch Uncuffed™ Murphy Eye ETT, preloaded with Slick Stylette® Device
- Preassembled and ready to use, saving time and space
- Cork-style stopper helps maintain stylet position

[View Details](#)



RUSCH® SAFETY CLEAR® UNCUFFED ETT

[View Details](#)



RUSCH® SAFETY CLEAR® PLUS CUFFED ETT

[View Details](#)



SHERIDAN® EZ-ENDO® ENDOTRACHEAL TUBE

The Sheridan® EZ-Endo® Endotracheal Tube features:

- Preassembled and ready to use, saving time and space
- Single-use, sterile
- Not made with natural rubber latex

[View Details](#)

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.