

Home // Airway Management // Endotracheal Tubes // Stylets

## STYLETS

Teleflex provides several Rusch® and Sheridan® Stylets to facilitate your intubation procedure.



### RUSCH® FLEXI-SLIP™ STYLET

The Rusch® Flexi-Slip™ Stylet features:


- Plastic coated malleable metal core
- Soft, retracted metal tip
- Coated for easy insertion and withdrawal

[View Details](#)

### RUSCH® SLICK STYLETTE® DEVICE

The Rusch® Slick Stylette® Device features:


- Made of malleable nylon-coated wire that holds its shape to facilitate intubation
- Specially formulated coating does not crack or peel, while allowing for smooth tube insertion

[View Details](#)

### SHERIDAN® SHER-I-SLIP® STYLET

The Sheridan® Sher-I-Slip® Stylet features:

- Designed to retain desired curve during oral intubation
- Material that reduces stylet-to-tube friction

[View Details](#)

### SHERIDAN® SOFT-TIP™ STYLET

The Sheridan® Soft-Tip™ Stylet features:

- Soft tip created by the retracted metal core
- Slippery outer finish

[View Details](#)

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.