


CATALOG # : 1844

RESPIRATORY

Description:

Softech® Bi-Flo ETCO<sub>2</sub> Sampling Cannula, Adult, 25/CS (Individually Packaged)



To order this product, please call:

800-387-9699

Canada

Lat

Am

- Details:
- Adult Cannula with 7 ft oxygen supply tubing and 7 ft CO<sub>2</sub> sampling line with a male luer connector

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.