

## CONTINUOUS NERVE BLOCK CATHETERS

Teleflex is dedicated to providing a portfolio of continuous peripheral nerve block (cPNB) solutions that promote positive patient outcomes and maximize clinician efficiency. These Arrow® Products include cPNB catheters, needles and trays in a variety of configurations.



### FLEXBLOCK® CPNB CATHETER

The Arrow® FlexBlock® Continuous Peripheral Nerve Block (cPNB) Catheter is an echogenic, flexible, kink-resistant catheter designed to help clinicians better control patient pain and improve overall block success.

[View Details](#)

### STIMUCATH® CPNB CATHETER

The Arrow® StimuCath® Continuous Peripheral Nerve Block (cPNB) Catheter allows clinicians to objectively confirm accurate catheter tip placement through nerve stimulation, which can help improve secondary block success. The wire-reinforced catheter body and coil tip are visible under ultrasound, providing the clinician with additional feedback on the location of the catheter.

[View Details](#)

### NEEDLES FOR CPNB CATHETERS

Teleflex offers Arrow® cPNB Needles in a range of lengths and sizes for use with the Arrow® StimuCath® cPNB Catheter.

[View Details](#)