

## EPIDURALS

### SINGLE SHOT EPIDURAL KITS AND SETS

#### Description:

Teleflex provides a range of options for single shot epidurals, in kit and set configurations, including a bullet-tip epidural needle designed to deflect when it contacts the dura. Insertion is aided by an introducer needle.



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

#### Intermediate Sets

#### Kits

© 2021 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.