

CONTINUOUS NERVE BLOCK CATHETERS**STIMUCATH® CPNB CATHETER****Description:**

The Arrow® StimuCath® Continuous Peripheral Nerve Block (cPNB) Catheter allows clinicians to objectively confirm accurate catheter tip placement through nerve stimulation, which can help improve secondary block success. The wire-reinforced catheter body and coil tip are visible under ultrasound, providing the clinician with additional feedback on the location of the catheter.

All StimuCath® Trays include:

- StimuCath® Continuous Peripheral Nerve Block Catheter
- Tuohy needle: specially designed to allow clinicians to retract the catheter, if necessary, back into the needle without withdrawing the needle and resticking the patient (see IFU for complete instructions).
- Stimulating SnapLock™ Catheter-Syringe Adapter
- StatLock® Stabilization Device

Safety Kits Include:

- 3mL ChloraPrep® with Hi-Lite Orange® Tint
- SnapTear® Drape with 5" fenestration
- Safety Injection Needles
- Procedural medications
- Procedural components and dressings

Kits Include:

- Procedural medications
- Procedural components and dressings

Sets Include:

- Basic procedural components



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Safety Kits

Sets

© 2021 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.