

Home // Pain Management // Peripheral Nerve Block Needles // UltraQuik® Tuohy PNB Needle (echogenic)

PERIPHERAL NERVE BLOCK NEEDLES

ULTRAQUIK® TUOHY PNB NEEDLE (ECHOGENIC)

Description:

The Arrow® UltraQuik® Tuohy Echogenic PNB Needle is for clinicians who prefer to only use ultrasound guidance when performing single-shot peripheral nerve blocks. The UltraQuik® Tuohy Needle can help increase overall block success by enhancing needle visibility under ultrasound, and needle control and maneuverability during needle positioning.

- Five grooved rings at the distal end of the needle cannula to help clinicians identify the needle tip under ultrasound
- Designed with tuohy tip for optimal balance between tactile feel and penetration
- Uninsulated needle shaft for ultrasound-only needle guidance
- Centimeter markings on needle shaft for visual confirmation of insertion depth
- Ergonomic needle hub for precision grip
- Extended injection tubing length helps reduce the risk of contamination to the sterile field



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

UltraQuik® Tuohy PNB Needle (echogenic)

© 2022 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.