

CATALOG # : AK-05000

THERACATH® EPIDURAL KITS AND SETS**Description:**

Arrow® Epidural Catheterization Intermediate Set with TheraCath® Catheter with Stylet, 10/CS



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Contents:

- 1: Epidural Catheter: 19 Ga. TheraCath® with Stylet, SnapLock™ Adapter
- 1: Epidural Needle: 17 Ga. x 3-7/8" (9.84 cm) TW with cm markings
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 5 cc Luer-Lock Loss of Resistance, Glass
- 1: Syringe: 20 cc Luer-Lock
- 1: Syringe: 3 cc Luer-Lock
- 1: Medicine Cup
- 1: Drape: 24" x 36"
- 1: Drape: 24" x 30" with 4" fenestration, with adhesive
- 1: CSR Wrap
- 1: Filter: 5 Micron Straw
- 4: Gauze Pads: 4" x 4"
- 2: Prep Sponge Swabs
- 1: Wing Clip: Hub
- 1: Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.