

CATALOG # : ASA-27090-SB

**SUREBLOCK® SPINAL KITS AND SETS****Description:**

Arrow® Spinal Anesthesia Kit with 27 Ga. x 3-1/2" Sprotte® Spinal Needle, 5/CS



United States

To order this product, please call:

**866-246-6990**

Lat

Am

**Contents:**

- 1: Sprotte® Spinal Needle: 27 Ga. x 3-1/2" (90 mm) with 22 Ga. x 30 mm Introducer
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip, Glass Barrel
- 1: 1 mL Ampule 0.1% Epinephrine Solution
- 1: 2 mL Ampule 0.75% Bupivacaine and 8.25% Dextrose Solution
- 1: 5 mL Ampule 1% Lidocaine HCl Solution
- 1: Foil Pack 10% Povidone-Iodine Solution
- 1: SharpsAway® Disposal Cup
- 1: CSR Wrap
- 1: Clear Fenestrated Drape
- 1: Towel
- 1: Particulate Filter: 5 Micron Straw
- 1: Medication Label Sheet
- 4: Gauze Pads: 4" x 4"
- 3: Prep Sponge Swabs
- 1: Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.