


CATALOG # : FB-19604-SST

FLEXBLOCK® CPNB CATHETER

Description:

Arrow® FlexBlock® 19 Ga. x 60 cm Echogenic cPNB Catheter Set with 17 Ga. x 4 cm Stimulating Echogenic Tuohy, 10/CS



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Contents:

- 1: Continuous Peripheral Nerve Block Catheter: Echogenic 19 Ga. (1.1 mm) x 60 cm FlexBlock® with Stylet and Threading Assist Device
- 1: Stimulating Plexus Block Needle: Echogenic 17 Ga. x 1-1/2" (4 cm) TWS, with cm markings
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Skin Protectant Prep Pad
- 1: Catheter Syringe Adapter: SnapLock™
- 1: Filter: 0.2 Micron Flat
- 1: Dressing: StatLock® Catheter Stabilization Device 2-1/2" x 3-3/4" (6.7 x 9.5 cm) with 6" (15 cm) Extension Tubing
- 1: Tubing: 18" (46 cm) Extension

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.