

CATALOG # : FB-19605-SNS

**FLEXBLOCK® CPNB CATHETER****Description:**

Arrow® FlexBlock® 19 Ga. x 60 cm Echogenic cPNB Catheter Set with 17 Ga. x 5 cm Echogenic Tuohy, 10/CS



United States

To order this product, please call:

**866-246-6990**

Canada

To order this product, please call:

**800-387-9699**

Lat

Am

**Contents:**

- 1: Continuous Peripheral Nerve Block Catheter: Echogenic 19 Ga. (1.1 mm) x 60 cm FlexBlock® with Stylet and Threading Assist Device
- 1: Plexus Block Needle: Echogenic 17 Ga. x 2" (5 cm) TWS, with cm markings
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Skin Protectant Prep Pad
- 1: Catheter Syringe Adapter: SnapLock™
- 1: Filter: 0.2 Micron Flat
- 1: Dressing: StatLock® Catheter Stabilization Device 2-1/2" x 3-3/4" (6.7 x 9.5 cm) with 6" (15 cm) Extension Tubing
- 1: Tubing: 18" (46 cm) Extension

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.