


CATALOG # : FB-19608-SST

FLEXBLOCK® CPNB CATHETER

Description:

Arrow® FlexBlock® 19 Ga. x 60 cm Echogenic cPNB Catheter Set with 17 Ga. x 8 cm Stimulating Echogenic Tuohy, 10/CS



United States

To order this product, please call:
866-246-6990

Canada

To order this product, please call:
800-387-9699

Lat
Am

- Contents:
- 1: Continuous Peripheral Nerve Block Catheter: Echogenic 19 Ga. (1.1 mm) x 60 cm FlexBlock® with Stylet and Threading Assist Device
 - 1: Stimulating Plexus Block Needle: Echogenic 17 Ga x 3" (8 cm) TWS, with cm markings
 - 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
 - 1: Skin Protectant Prep Pad
 - 1: Catheter Syringe Adapter: SnapLock™
 - 1: Filter: 0.2 Micron Flat
 - 1: Dressing: StatLock® Catheter Stabilization Device 2-1/2" x 3-3/4" (6.7 x 9.5 cm) with 6" (15 cm) Extension Tubing
 - 1: Tubing: 18" (46 cm) Extension