

CATALOG # : FB-19611-K

FLEXBLOCK® CPNB CATHETER

Description:

Arrow® FlexBlock® 19 Ga. x 60 cm Echogenic cPNB Catheter Kit with 17 Ga. x 11 cm Stimulating Echogenic Tuohy, 5/CS



United States

To order this product, please call:

866-246-6990

Lat

Am

Contents:

- 1: Continuous Peripheral Nerve Block Catheter: Echogenic 19 Ga. (1.1 mm) x 60 cm FlexBlock® with Stylet and Threading Assist Device
- 1: Stimulating Plexus Block Needle: Echogenic 17 Ga. x 4-3/8" (11 cm) TWS with cm markings
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 20 mL Luer-Lock
- 1: Syringe: 3 mL Luer-Lock
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep® One-Step Solution with Hi-Lite Orange® Tint
- 1: 5 mL Ampule 1% Lidocaine HCl Solution
- 1: 5 mL Ampule 1.5% Lidocaine HCl with Epinephrine 1:200,000 Solution
- 1: Skin Protectant Prep Pad
- 1: Swabstick Foil Pack Skin Protectant
- 1: Catheter Syringe Adapter: SnapLock™
- 1: SharpsAway® II Locking Disposal Cup
- 1: CSR Wrap
- 1: Drape: SnapTear® Clear 30" x 26" with 5" fenestration
- 1: Towel
- 1: Filter: 0.2 Micron Flat
- 1: Particulate Filter: 5 Micron Straw
- 1: Medication Label: 1% Lidocaine
- 1: Medication Label: 1.5% Lidocaine with Epinephrine 1:200,000
- 1: Dressing: StatLock® Catheter Stabilization Device 6.7 cm x 9.5 cm with 15 cm Extension Tubing
- 2: Dressings: Tegaderm® 10 cm x 12 cm

- **8:** Gauze Pads: 4" x 4"
- **1:** Tubing: 18" (46 cm) Extension

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.