

CATALOG # : JC-05400-E

FLEXTIP PLUS® EPIDURAL KITS AND SETS

Description:

Arrow® Epidural Catheterization Set with FlexTip Plus® Open Tip Catheter, 20/CS



United States

To order this product, please call:  
**866-246-6990**

Canada

To order this product, please call:  
**800-387-9699**

Lat

Am

- Contents:
- 1: Epidural Catheter: 19 Ga. x 90 cm FlexTip Plus® (single open end hole) with SnapLock™ Adapter, Threading Assist Device
  - 1: Epidural Needle: 17 Ga. x 3-7/8" (9.84 cm) TW with cm markings
  - 1: Syringe: 10 cc Luer-Slip Loss of Resistance
  - 1: Wing Clip: Hub
  - 1: Filter: 0.2 Micron Flat

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.