

CATALOG # : MH-05001

THERACATH® EPIDURAL KITS AND SETS**Description:**

Arrow® Epidural Catheterization Kit with TheraCath® Catheter with Stylet, 10/CS



United States

To order this product, please call:

866-246-6990

Lat

Am

Contents:

- 1: Epidural Catheter: 19 Ga. TheraCath® with Stylet, SnapLock™ Adapter
- 1: Epidural Needle: 17 Ga. x 3-7/8" (9.84 cm) TW with cm markings
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 20 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip Loss of Resistance, Glass
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: 10 mL Ampule 0.9% Saline Solution
- 1: 5 mL Ampule HCl, 1% Lidocaine Solution
- 1: 5 mL Ampule HCl, 1.5% Lidocaine w/Epinephrine 1:200,000 Solution
- 1: Wing Clip: Hub
- 1: SharpsAway® II Locking Disposal Cup
- 1: CSR Wrap
- 1: Drape: 24" x 36"
- 1: Drape: 24" x 40" with 2-5/8" fenestration, with adhesive
- 1: Filter: 0.2 Micron Flat
- 1: Filter: 5 Micron Straw
- 1: Epidural Pad
- 1: Medication Label: 0.9% Saline
- 1: Medication Label: 1% Lidocaine
- 1: Medication Label: 1.5% Lidocaine with Epinephrine 1:200,000
- 8: Gauze Pads: 4" x 4"
- 3: Prep Sponge Swabs
- 1: Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.