

CATALOG # : MP-17019-TIG

FLEXTIP PLUS® EPIDURAL KITS AND SETS

Description:

Arrow® Epidural Catheterization Intermediate Set with FlexTip Plus®, Closed Tip Multi-Port Catheter, 5/CS



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat

Am

Contents:

- 1: Epidural Catheter: 19 Ga. x 90 cm FlexTip Plus®, Closed Tip, Multi-Port (4 lateral holes) with Threading Assist Device
- 1: Epidural Needle: Tuohy 17 Ga. x 3-1/2" (8.89 cm) TW with cm markings
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 20 mL Luer-Lock
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Lock Loss of Resistance, Glass
- 1: Catheter Syringe Adapter: SnapLock™ Adapter
- 1: SharpsAway® II Locking Disposal Cup
- 1: CSR Wrap
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 30" x 24" (76 cm x 61 cm) with 4" (10 cm) fenestration and adhesive
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Filter: 0.2 Micron Flat
- 1: Particulate Filter: 5 Micron Straw
- 1: Epidural Pad
- 1: Dressing: Tegaderm® 10 cm x 12 cm
- 5: Gauze Pads: 4" x 4" (10 cm x 10 cm)
- 3: Prep Sponge Swabs
- 1: Medicine Cup
- 1: Tray: Prep

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.