

CATALOG # : MT2060LYK5-CP

AUTOFUSER® PAIN PUMP FOR INCISIONAL INFUSIONS

Description:

Arrow® AutoFuser® Pain Pump for Incisional Infusions: 275 x 2 mL/hr, 2 mL-60 min Refill, with YK5 Catheter Kit, 5/BX



United States

To order this product, please call:  
**866-246-6990**

Lat  
Am

Contents:

- **Catheter Kit includes:**
- **2:** 5" (12.8 cm) Saturation Catheters
- **2:** 4 Fr x 4" (10.2 cm) Peel Away Introducer Needles
- **1:** Bifurcated Set
- **2:** Occlusive Dressings
- **3:** Adhesive Strip Cards
- **2:** Prime/Flush Syringes
- **2:** Luer Caps
- **1:** Prep Tray
- **1:** Filling Syringe