

Home // Pain Management // Airways // ETCO₂ // Softech® Bi-Flo ETCO₂ Sampling Cannula // RHET43U

CATALOG # : RHET43U

RESPIRATORY

Description:

Softech® Bi-Flo ETCO₂ Sampling Cannula, Adult, 25/CS (Individually Packaged)



United States

To order this product, please call:

866-246-6990

Contents:

- Adult Cannula with 7 ft oxygen supply tubing and 3 in. CO₂ sampling line with a female luer connector

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.