

CATALOG # : SL-05000

**THERACATH® EPIDURAL KITS AND SETS****Description:**

Arrow® Epidural Catheterization Intermediate Set with TheraCath® Catheter with Stylet, 10/CS



United States

To order this product, please call:

**866-246-6990**

Lat

Am

**Contents:**

- 1: Epidural Catheter: 19 Ga. TheraCath® with Stylet, SnapLock™ Adapter
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 16 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Epidural Needle: 17 Ga. x 3-7/8" (9.84 cm) TW with cm markings
- 1: Syringe: 5 cc Loss of Resistance, Glass
- 1: Syringe: 20 cc Luer-Slip
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 cc Luer-Lock Syringe
- 1: Foil Pack: 10% Povidone-Iodine Solution
- 1: Medicine Cup
- 1: CSR Wrap
- 1: Drape: 24" x 36" with 4" fenestration, with adhesive
- 1: Filter: 5 Micron Straw
- 1: Epidural Pad
- 8: Gauze Pads: 4" x 4"
- 3: Prep Sponge Swabs
- 1: Wing Clip: Hub
- 1: Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.