

CATALOG # : SM-03001

SINGLE SHOT EPIDURAL KITS AND SETS

Description:

Arrow® Single Shot Epidural Kit with Tuohy Needle, 10/CS



United States

To order this product, please call:
866-246-6990

Lat
Am

Contents:

- 1: Epidural Needle: 18 Ga. x 3-1/2" (8.89 cm) TW with cm markings
- 1: Filter Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 5 cc Luer-Lock Loss of Resistance, Glass
- 1: Syringe: 10 cc Luer-Lock
- 1: Syringe: 3 cc Luer-Lock
- 1: Syringe: 3 cc Luer-Slip
- 1: Foil Pack: 10% Povidone-Iodine Solution
- 1: 10 ml Ampule 0.9% Saline Solution
- 2: 5 ml Ampules HCl, 1% Lidocaine Solution
- 1: Adapter Tubing: Female Luer-Lock/Male Luer-Slip
- 1: Towel: 13" x 18"
- 1: Drape: 24" x 30" with 4" fenestration, with adhesive
- 1: CSR Wrap
- 5: Gauze Pads: 4" x 4"
- 3: Prep Sponge Swabs
- 1: Tray: Prep