

CATALOG # : SM-03001

SINGLE SHOT EPIDURAL KITS AND SETS**Description:**

Arrow® Single Shot Epidural Kit with Tuohy Needle, 10/CS



United States

To order this product, please call:

866-246-6990

Lat

Am

Contents:

- 1: Epidural Needle: 18 Ga. x 3-1/2" (8.89 cm) TW with cm markings
- 1: Filter Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 5 cc Luer-Lock Loss of Resistance, Glass
- 1: Syringe: 10 cc Luer-Lock
- 1: Syringe: 3 cc Luer-Lock
- 1: Syringe: 3 cc Luer-Slip
- 1: Foil Pack: 10% Povidone-Iodine Solution
- 1: 10 ml Ampule 0.9% Saline Solution
- 2: 5 ml Ampules HCl, 1% Lidocaine Solution
- 1: Adapter Tubing: Female Luer-Lock/Male Luer-Slip
- 1: Towel: 13" x 18"
- 1: Drape: 24" x 30" with 4" fenestration, with adhesive
- 1: CSR Wrap
- 5: Gauze Pads: 4" x 4"
- 3: Prep Sponge Swabs
- 1: Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.