

CATALOG # : WJ-05401

FLEXTIP PLUS® EPIDURAL KITS AND SETS

Description:

Arrow® Epidural Catheterization Kit with FlexTip Plus® Open Tip Catheter, 10/CS



United States

To order this product, please call:

866-246-6990

Lat

Am

Contents:

- 1: Epidural Catheter: 19 Ga. x 90 cm FlexTip Plus® (single open end hole) with Threading Assist Device
- 1: Epidural Needle: Hustead 17 Ga. x 3-7/8" (9.84 cm) TW with cm markings
- 1: Injection Needle: 16 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 20 mL Luer-Slip
- 1: Syringe: 5 mL Luer-Slip Loss of Resistance, Glass
- 1: Injection Needle: 25 Ga. x 5/8" (1.59 cm) and 3 mL Luer-Lock Syringe
- 1: 10 mL Ampule 0.9% Saline Solution
- 1: 5 mL Ampule 1% Lidocaine HCl Solution
- 1: 5 mL Ampule 1.5% Lidocaine HCl with Epinephrine 1:200,000 Solution
- 1: Foil Pack 10% Povidone-Iodine Solution
- 1: Catheter Syringe Adapter: SnapLock™ Adapter
- 1: Wing Clip: Hub
- 1: SharpsAway® II Locking Disposal Cup
- 1: Clear Fenestrated Drape with adhesive
- 1: Towel: 13" x 18"
- 1: CSR Wrap
- 1: Filter: 0.2 Micron Flat
- 1: Particulate Filter: 5 Micron Straw
- 1: Epidural Pad
- 1: Medication Label: 0.9% Saline
- 1: Medication Label: 1% Lidocaine
- 1: Medication Label: 1.5% Lidocaine with Epinephrine 1:200,000

- **8:** Gauze Pads: 4" x 4"
- **3:** Prep Sponge Swabs
- **1:** Medicine Cup
- **1:** Tray: Prep

© 2024 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions. Contact customer service to confirm availability in your region.