

MAXTEC BLENDERS AND ACCESSORIES

BLENDERS

Description:

The Maxtec® MaxVenturi® Blender Setup delivers up to 60 LPM and provides for quick, comfortable, and effective patient use. It only requires piped oxygen and uses the Venturi principle to mix ambient air and oxygen without the need for piped medical grade air.

The MicroMax™ Blenders require oxygen and medical air connection and features delivery of accurate FiO₂ mixtures from two outlet ports. The High Flow version has a primary outlet flow range of 15-120 LPM. The Low Flow version has a primary outlet flow range of 3-30 LPM. When paired with a Maxtec® Blender Buddy, the MircoMax® provides a total mixing package, complete with blender, two flow meters, and an oxygen analyzer.*



United States

To order this product, please call:

866-246-6990

Canada

To order this product, please call:

800-387-9699

Lat Am

Catalog #	Description	Qty/Case
R211P03-036	MaxVenturi® Blender	1
R203P14	MicroMax™ High Flow	1
R203P13	MicroMax™ Low Flow	1

* The Hudson RCI® Neptune® Heated Humidifier is indicated for use at flows of 1-60 LPM when used with the Comfort Flo® System and Comfort Flo® Plus Cannula. Refer to the Hudson RCI® Neptune® Heated Humidifier user manual for more information.

©2019 Teleflex Incorporated. All rights reserved.

The listed Teleflex® item(s) are not in every case completely identical to the item illustrated but are the closest match available. Not all products are available in all regions; contact customer service to confirm availability in your region.